

GGM NP SERIES INFANT BUBBLE CPAP SYSTEM




USE IN NEONATE CPAP AND BUBBLE CPAP SYSTEM

GGM NP SERIES INFANT BUBBLE CPAP SYSTEM, COMPLETE SYSTEM INCLUDING SERVO CONTROL HUMIDIFIER, FULL RANGES OF NASAL PRONG, BONNET, CPAP GENERATOR, HEATER WIRE CIRCUIT, PRESSURE RELEASE SAFETY VALVE.

COMPLETE SYSTEM FOR APPLICATION OF CONTINUOUS POSITIVE AIRWAY PRESSURE TO NEONATES AND INFANTS BY NASAL PRONGS (NCPAP) ADMINISTERED WITH A HEATER WIRE CIRCUIT USED IN CONJUNCTION WITH A CONTINUOUS FLOW SOURCE INFANT VENTILATOR OR A SUITABLY EQUIPPED MULTIPURPOSE VENTILATOR.

ORDER INFORMATION

PN	Description	PK
VN-BC01	Bubble CPAP System -1x Bubble CPAP Circuit w/chamber -1x Safety Valve Kit -1x CPAP Generator	20 sets/ctn
VN-BC02	Bubble CPAP System -1x VH-2600A Servo Control Humidifier w/Heater wire adaptor & Temp probe -3x Bubble CPAP Circuit w/chamber -3x Safety Valve Kit -3x CPAP Generator -3x Nasal prong Kit (All size:0~5) -3x Bonnet (All size:SML)	3 sets/ctn



PN	Description	PK
VN-BD01	Air/O2 Blender 2-18 LPM -1x Air pressure line -1x O2 pressure line -1x C-Clamp	8 sets/ctn
VN-BD02	Air/O2 Blender 5-50 LPM -1x Air pressure line -1x O2 pressure line -1x C-Clamp	8 sets/ctn

PN	Description	PK
VN-0001	Nasal prong kit #0 (7.5fr) <700g	10 sets/box
VN-0002	Nasal prong kit #1 (10fr) 700~1250g	10 sets/box
VN-0003	Nasal prong kit #2 (12.5fr) 1250~2000g	10 sets/box
VN-0004	Nasal prong kit #3 (15fr) 2000~3000g	10 sets/box
VN-0005	Nasal prong kit #4 (16fr) <3000g	10 sets/box
VN-0006	Nasal prong kit #5 (20fr) 1~2 years	10 sets/box

PN	Description	PK
VN-HC01	Head cap / Bonnet #S	10 sets/box
VN-HC02	Head cap / Bonnet #M	10 sets/box
VN-HC03	Head cap / Bonnet #L	10 sets/box

PN	Description	PK
VN-BD03	Cart / Trolley with support arm with water bag arm	1 set/ctn



Great Group Medical Co., Ltd.

No.168,Xingong 2nd Rd.,Tianzhong Township,Changhua County,Taiwan

TEL:886-4-8758181 FAX:886-4-8756161

www.greatgroup.com.tw

info@greatgroup.com.tw